

B6.Cg-Tg(Cr2-cre)3Cgn/J Genotyping

1. Prepare DNA by standard isopropanol/ethanol extraction.
2. Bring final volume of each oligo to 100 pmol/ μ l (each) using ultrapure H₂O.

wt	KMO280 CTAGGCCACAGAATTGAAAGATCT KMO281 GTAGGTGGAAATTCTAGCATCATCC
mut	KMO278 GCGGTCTGGCAGTAAAACTATC KMO279 GTGAAACAGCATTGCTGCTCACTT

3. Set up a master mix (on ice) for each primer as follows:

12.5 μ L	Promega GoTaq Mastermix
.25 μ L	Primer Mix (=10pmol each primer)
19.8 μ L	H ₂ O



X No. of samples

4. Add 1 μ l template DNA to corresponding well in PCR plate.
5. Add 24 μ l master mix to each well in PCR plate
5. Seal using Microseal film.
6. Keep reactions on ice and load onto PCR block pre-heated to 94°.

PCR Program

94° 3 min.	} Repeat steps 2-4 for 35 cycles
94° :30 sec	
51.7° 1 min	
72° 1 min	
72° 2 min	
10° forever	

Expected band sizes:

WT: approx. 324bp
Tg: approx. 100bp